


Software Version : 6.3.1.0504
Sample Name : GHK-Cu brazowa fiolka
Instrument Name : DIONEX
Rack/Vial : 0/0
Sample Amount : 1,000000
Cycle : 1

Date : 2023-03-18 18:09:54
Data Acquisition Time : 2005-02-14 11:27:36
Channel : B
Operator : manager
Dilution Factor : 1,000000

Result File :
Sequence File : C:\PenExe\TcWS\Ver6.3.1\Examples\IC-106-S\ROK 2023\LUTY\14.02-peptydy\14.02-ghk.seq



REPORT

Peak #	Component Name	Time [min]	Amount [mg/fiolka]
4	GHK-Cu	3,678	23,79

Eko-Analiza M.Turek, M.Wysocki sp.j.
ul. Harnasiow 9; 43-300 Bielsko-Biala
eko_analiza@onet.pl